

“DAIHAN” **WiseMix® RK** Programmable Digital Rockers, Programmable Digital Feedback Control Smooth Start & Wave Motion, High Precision Digital Speed Control, 5~50 rpm, Tilt Angle $\pm 6 / \pm 9$



Programmable Digital Rocker (RK-1D)

- High Precision Digital Speed Control : 1rpm
- **Patented Jog-Shuttle Control System**
- CE Certified
- **“PL(Product Liability) Insurance”**
: One Hundred Million/Korean Won(₩) for Each of Human and an Object

Application & Features

- * High Precision Digital Speed Control : 1rpm
- * Programmable Operation Mode Supported : 10 Steps and 6 Programs
- * Digital Timer : Max. 99hr 59min
- * Digital LCD Display with Back-Light Function for Timer, RPM and Power Bar-Graph
- * Smooth Start and Change of the Stirring Speed
- * Over-Load Protection Functions for Motor
- * Storage Function for Set Values : RPM and Timer
- * Alarm Function : Error Status and Timer-end
- * Stainless Steel Platform Comes with Gray Nonslip Rubber Pad
- * Stackable Platform with Pad to Double Mixing Area (Optional)
- * Locking Mode supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)
- * **Patented Jog-Shuttle Control System**
- * CE Certified
- * **“PL(Product Liability) Insurance”**
: One Hundred Million/Korean Won(₩) for Each of Human and an Object

Patented Control System (No. 10-0554098)

1 or 2 Mode Buttons + Jog Shuttle Switch
Realize “Easy to Use”

Specification

Model		RK-1D	RK-2D
Dimensions (w×d×h)	Platform (mm)	290×200	300×300
	Overall (mm)	290×240×150	300×330×150
Angle (degree)		$\pm 6 / \pm 9$ (Default)	
Motor		35W Geared DC Motor	
Speed Range (rpm)		5 ~ 50	
Controller		Programmable Digital Feedback Controller with Jog-Shuttle Switch(Turn+Push)	
Program Function		10 Steps, 6 Programs	
Display		Digital LCD Display with Back-Light Function	
Timer		99hr 59min (Programmable)	
Safety Device		Over Load Protector for motor	
Others		CE Certified PL (Product Liability) Insurance Smooth Start and Change of the Stirring Speed for Stable Running Storage Function (RPM, Timer) Locking Mode (Jog-Shuttle Input Disabled) Alarm (Error Status and Timer-end)	
Packing Dimension (w×d×h) (mm)		450×350×250	
Gross Weight		7 kg	8 kg
Power Consumption		35 W	
Power Supply*		Free Voltage (AC 80V~240V, 50/60Hz)	

* Other Specifications are available upon Customer's Request.

Ordering Information

	Ea.	Price, ₩
DH.WRK01010 Digital Rocker, RK-1D , w290×d200mm Platform, with Rubber Pad		478,290
DH.WRK01020 Digital Rocker, RK-2D , w300×d300mm Platform, with Rubber Pad		562,190

Optional Accessories

	Ea.	Price, ₩
DH.WRK05011 Second Platform, RK5011 , for RK-1D , with 4 Pillars, w290×d200		42,920
DH.WRK05021 Second Platform, RK5021 , for RK-2D , with 4 Pillars, w300×d300		51,190