


Hotplate Magnetic Stirrer



MYP11-2



MYP11-2A



SH05-3G



X85-2



X85-2S

Technical Parameters :

| Model | MYP11-2 | MYP11-2A | SH05-3G | X85-2 | X85-2S |
|------------------------|--------------------------|---------------|-----------------|----------------|--------|
| Work Plate Size | φ145mm | | φ175mm | φ145mm | |
| Work Plate Material | Aluminum alloy | | Stainless steel | Aluminum alloy | |
| Max. Stirring Capacity | 10L | | | 5L | |
| Motor Type | Brushless DC motor | | | | |
| Motor Output | 15W | | | 10W | |
| Speed Range | 50-1500rpm | | 0-1500rpm | | |
| Digital Display | Speed, temperature, time | | | / | Speed |
| Control | Button | Knob + Button | Button | Knob | |
| Heating Output | 600W | | 800W | 600W | |
| Heating Temp. Range | RT~320°C | | RT~250°C | RT~300°C | |
| Temp. Accuracy | ±0.5°C | | | / | |
| Timing Range | 1~9999min | 1~250min | / | | |
| Power Supply | AC200V, 50/60Hz | | | | |
| External Size(W*D*H) | 280*160*90mm | | 330*190*130mm | 250*150*135mm | |
| Package Size(W*D*H) | 380*300*190mm | | 430*330*230mm | 350*290*240mm | |
| Gross Weight | 4.5kg | | 7kg | 6kg | |